1 Title: DELT-Hit: An end-to-end computational framework for DNA-encoded chemical library analysis

1.1 Authors

Adriano Martinelli, Alice Lessing, Garry Hoppeler, Andreas Glogger, Jörg Scheuermann

1.2 Abstract

DNA-encoded chemical libraries (DELs) have revolutionized drug discovery by enabling the simultaneous screening of millions to billions of small molecules through DNA-tag identification via high-throughput sequencing. As outlined in the comprehensive Nature Reviews Methods Primers on this technology (Satz et al., 2022), DELs are now employed by numerous pharmaceutical companies and academic laboratories worldwide. However, the computational analysis of DEL screening data remains a critical bottleneck, requiring sophisticated integration of genomics, cheminformatics, and statistical analysis workflows that are currently accessible only through proprietary or highly specialized software solutions. This protocol presents DELT-Hit (DNA-Encoded Library Technology Hit), a comprehensive open-source computational framework that makes DEL data analysis accessible through an intuitive command-line interface to both computational and experimental researchers. DELT-Hit is specifically designed to handle the scale and complexity of modern industrial DEL campaigns, supporting libraries containing hundreds of millions of compounds while maintaining computational efficiency and user accessibility. DELT-Hit offers a complete pipeline that converts raw FASTQ reads into machine learning ready chemical information through five connected modules: (1) adaptive sequence demultiplexing using optimized RNA-seq algorithms with DEL-specific error correction and flexible barcode handling, (2) automated chemical structure reconstruction from building block libraries using reaction SMARTS templates with support for both single and dual display architectures, (3) comprehensive molecular property calculation and descriptor generation using established cheminformatics libraries, (4) statistical analysis and hit ranking with multiple normalization strategies adapted from proven RNA-seq methodologies, and (5) integrated quality control and visualization tools specifically designed for DEL data interpretation. The modular architecture allows researchers to customize workflows while maintaining reproducibility through configuration files and standardized output formats. We demonstrate the protocol's effectiveness using representative single and dual display DEL screening datasets, showcasing the complete analysis pipeline from raw sequencing reads to ranked lists of chemical hits with computed chemical properties and representations for downstream machine learning tasks. The entire analysis, including quality control and visualization, can be completed within 4-6 hours on standard computational hardware for typical datasets, making it accessible to laboratories without specialized computing infrastructure. DELT-Hit addresses the critical computational gap identified in the DEL field and provides the standardization

necessary for reproducible analysis across research groups. The protocol is accompanied by comprehensive documentation, tutorial datasets with both single and dual display examples, ensuring broad adoption and consistent implementation across the growing DEL community.

1.3 Key Points

- Industrial-scale capabilities: DELT-Hit is designed to handle the computational demands of modern pharmaceutical DEL campaigns, efficiently processing libraries containing hundreds of millions of compounds while maintaining user-friendly operation
- Comprehensive dual architecture support: The framework provides native support for both single and dual display DEL architectures, addressing the full spectrum of current library designs used in industry and academia
- Validated algorithms: Integrates proven bioinformatics tools (Cutadapt for sequence processing, edgeR for statistical analysis) with specialized DEL-specific optimizations, error handling, and quality control metrics developed through extensive validation studies
- Flexible and robust design: Modular architecture accommodates diverse library formats, custom reaction templates, and building block definitions
- Research-grade quality assurance: Built-in quality control metrics, automated validation checks, and standardized reporting ensure reliable results and facilitate systematic troubleshooting across different experimental conditions
- Machine learning ecosystem integration: Generates standardized, analysis-ready datasets fully compatible with downstream machine learning workflows for advanced hit prediction, structure-activity relationship analysis, and virtual screening applications

1.4 Technical Overview

DELT-Hit is implemented as a Python package organized into five modules:

1.5 Core Analysis Modules:

- init: Project initialization and configuration management
- demultiplex: Sequence processing and demultiplexing with adaptive error correction
 - qc: produce quality control plots
 - report: produce report with statistics about mapped sequences
- library: Chemical structure reconstruction and molecular property calculation
 - enumerate: construct smiles from reaction steps
 - properties: compute and visualize distribution of chemical properties
 - represent: compute SMILES representations for downstream tasks

• dashboard: interactive data exploration and visualization

• analyse: Statistical analysis and hit ranking

1.6 Introduction

DNA-encoded chemical libraries (DECLs) have emerged as a powerful technology in drug discovery, enabling the synthesis and screening of vast chemical spaces that would be impractical to explore using traditional approaches. In DECL technology, each chemical compound is covalently linked to a unique DNA barcode, allowing millions to billions of compounds to be screened simultaneously against biological targets through DNA sequencing-based identification. The computational analysis of DECL screening data presents unique challenges that require specialized approaches: accurate demultiplexing of complex DNA barcode combinations from sequencing reads, reconstruction of chemical structures from building block combinations, statistical analysis of enrichment patterns, and integration with cheminformatics workflows for hit optimization. Current computational tools for DECL analysis often focus on individual workflow components rather than providing comprehensive solutions, require significant programming expertise, or lack flexibility for diverse library architectures. Most existing approaches do not integrate well with standard bioinformatics pipelines or provide adequate quality control mechanisms.

1.6.1 Development of the protocol

DELT-Core addresses these limitations through a unified, modular framework built around several key design principles: leveraging established bioinformatics tools where appropriate, providing flexible configuration for diverse library designs, implementing comprehensive quality control, and maintaining accessibility for users with varying computational backgrounds. The framework consists of integrated modules: demultiplexing using adapted Cutadapt workflows, chemical structure reconstruction with RDKit, statistical analysis with edgeR, molecular property calculation, and interactive visualization dashboards. This modular architecture enables users to execute complete workflows or use individual components as needed. The protocol has been successfully applied to analyze diverse DECL architectures including multi-cycle libraries, hybridized libraries combining independent synthetic routes, and large-scale screens with millions of compounds. The framework's flexibility allows adaptation to novel experimental protocols with minimal code modification.

1.6.2 Comparison with other methods

A number of academic and commercial solutions address parts of the DEL informatics workflow, but few provide an end-to-end, openly available pipeline. DELT-Hit is most comparable in scope to recent academic offerings such as DELi (UNC), while also integrating component tools that specialize in key problem areas (e.g., demultiplexing with Cutadapt,

enrichment modeling with edgeR). Below we contrast DELT-Hit against representative methods across availability, scope/functionality, and performance/scalability.

1.6.3 Applications of the method

TODO: Add selected papers from group.

1.6.4 Limitations

- Computational requirements scale significantly with library size and sequencing depth
- Complex library architectures with non-standard reaction schemes may require customization
- Demultiplexing assumes independence of errors across barcode positions
- Chemical structure reconstruction depends on accurate reaction SMARTS definitions
- Very large datasets may require high-memory computing resources or distributed processing

1.6.5 Overview of the procedure

The protocol consists of five main stages executed through a collaborative command-line interface: (i) project initialization and library definition, (ii) chemical structure enumeration and property calculation, (iii) sequence demultiplexing and quality control, (iv) statistical analysis and hit identification, and (v) data visualization and interpretation. The workflow is designed to accommodate both computational biologists familiar with bioinformatics pipelines and medicinal chemists focused on chemical structure analysis. Each stage includes comprehensive quality control and generates standardized outputs compatible with downstream analysis tools.

1.7 Experimental design

1.7.1 Input requirements

The framework processes three primary input types: (1) library definition files containing building block structures, reaction SMARTS, and constant sequences; (2) experimental metadata specifying selection conditions and sample identifiers; and (3) raw FASTQ files from high-throughput sequencing platforms.

1.7.2 Library architecture considerations

The framework supports diverse library architectures:

- Single-display libraries
- Dual-display libraries

1.7.3 Library chemistry

• The reactions performed represented as a reaction graph. In theory, arbitrary reaction steps can be encoded and performed by the tool

1.7.4 Quality control parameters

Key quality metrics monitored throughout the workflow:

- Demultiplexing efficiency and barcode recovery rates
- Sequencing quality scores and adapter matching statistics
- Chemical structure validation and reaction success rates
- Statistical significance of enrichment patterns

1.8 Materials

1.8.1 Hardware:

- Minimum: 8 GB RAM, 8 CPU cores, 30 GB storage
- Recommended: 16 GB RAM, 32 CPU cores, 30 GB storage

1.8.2 Operating system:

• Linux (Ubuntu 22.04+), macOS (12.0+), or Windows

1.8.3 Software

- Python 3.10 or higher with conda package manager or virtual environment
- cutadapt (4.9+): Sequence adapter trimming and demultiplexing
- rdkit (2024.3+): Chemical structure processing and property calculation
- pandas (2.2+): Data manipulation and analysis
- matplotlib/seaborn: Data visualization and plotting
- edgeR (via R): Statistical analysis of count data

1.8.4 File Preparation

- 1. Sequencing file in formats compatible with cutadapt (<TODO: list formats />)
- 2. Configuration file

1.8.5 Setup

1.8.5.1 1. Conda We recommend using the Miniconda package manager to create an isolated environment for this project. This ensures that all dependencies are managed correctly.

- Download and install Miniconda for your operating system.
- After installation, you should be able to use the conda command in your terminal. bash conda create -n del python=3.11 -y conda activate del # Always activate this environment (`conda activate del`) before using `delt-core`. pip install git+https://github.com/DELTechnology/delt-core.git delt-cli --help # You should see a list of available commands.

1.8.5.2 2. R Environment Some analysis features in delt-core (like enrichment analysis with edgeR) depend on R.

- Install R: Download and install R from the Comprehensive R Archive Network (CRAN).
- Install R Packages: Once R is installed, open an R console and run the following commands to install the required packages: "'R # Install tidyverse and GGally from CRAN install.packages(c("tidyverse", "GGally"))

```
\# Install BiocManager if (!require("BiocManager", quietly = TRUE)) install.packages("BiocManager")
```

Install edgeR and limma from Bioconductor BiocManager::install(c("edgeR", "limma")) "'

1.9 Procedures

1. Create configuration file

Timing: 30 minutes - 8h

• Experiment section

variable	value
name	test-1
$fastq_path$	\sim /data/DECLT-DB/fastq_files/368061_1-
	$241105_AG_BZ_NC_pool1_NF_S3_R1_001.fastq.gz$
$save_dir$	\sim /data/DECLT-DB/experiments
num_cores	10

• Selections section

name opera	torlate	target	group	beads info	blocki	nguffer	protoco\$0	S1
AG24_A.	26-	No	no_pro	t Diy mabeads	Biotin	PBS-	DECL_A50Y	ACAGCI
Gloge	r Sep-	Protein		SA C1		Τ	_	
	24							
AG24_ A 2.	26-	No	no_pro	t Dy nabeads	Biotin	PBS-	DECL_A509	AGCACT
Gloge	r Sep-	Protein		SA C1		Τ		
	24							
$AG24$ _ A 3.	26-	No	no_pro	t Dy nabeads	Biotin	PBS-	DECL_A50	AT CG CT
Gloge	r Sep-	Protein		SA C1		Τ		
	24							
AG24_A0	26-	hCAII-	protein	Dynabeads	Biotin		DECL_A50	GACCC
Gloge	r Sep-	Avi-His		SA C1		Τ		
	24	(biot)						
AG24_A1	26-	hCAII-	protein	Dynabeads	Biotin		DECL_A50	GCCACI
Gloge	r Sep-	Avi-His		SA C1		Τ		
1001 40	24	(biot)		D 1 1	D	DDC		
AG24_A12	26-	hCAII-	protein	Dynabeads	Biotin		DECL_A50	MA G GC'I
Gloge	r Sep-	Avi-His		SA C1		Τ		
A C 0.4 M 0	24 26-	(biot) USP-13	:	D	D: -4:	DDC		T (1000
AG24_A19	r Sep-	(biot.)	naive	Dynabeads SA C1	DIOUII	Т	DECL_AGO	M COMPONE
Gioge	24	(1000.)		SA CI		1		
AG24_ A 0	26-	USP-13	naive	Dynabeads	Riotin	PRS_	DECL_A50	TCMCT
	r Sep-	(biot.)	naive	SA C1	Diomi	т Б Б-	DECL_AW	
Gloge	24	(5100.)		D11 O1		*		
AG24 A21	26-	USP-13	naive	Dynabeads	Biotin	PBS-	DECL_A503	ACTAC7
_	r Sep-	(biot.)		SA C1		T	_ = <u></u> 2G	
0-	24	()		-				

• Structure section

name	type	max_error_rate	indels
S0	selection	0	FALSE
C0	constant	0	FALSE
B0	building_block	0	FALSE
C1	constant	0	FALSE
B1	building_block	0	FALSE
C2	constant	0	FALSE
S1	selection	0	FALSE

name type max_error_rate	indels
--------------------------	--------

• Compounds section

name	smiles
scaffold_1	Ic1ccc(CC(N=[N+]=[N-])C(O)=O)cc1
scaffold_2	$[N-]{=}[N+]{=}NC(C(O){=}O)Cc1cc(I)ccc1$

• Reactions section

name	smirks
ABF	[CX3:1](=[0:2])[OX2;H1].[N;H2:4]>>[CX3:1](=[0:2])[N;H:4]
SR	[#6:1][\$([NX2-][NX2+]#[NX1]),\$([NX2]=[NX2+]=[NX1-])]>>[#6:1][N;H2]
CuAAC	$ [\texttt{CX2:1}] \# [\texttt{CX2:2}] \cdot [\texttt{N:3}] = [\texttt{N+:4}] = [\texttt{N-:5}] >> [\texttt{C:1}] \ 1 = [\texttt{C:2}] \ [\texttt{N:3}] \ [\texttt{N:4}] = [\texttt{N:5}] \ 1 $
Suz	[cX3:1][I].[#6:2][BX3]>>[cX3:1][#6:2]

• Constant section

name	codon
C0	GGAGCTTCTGAATTCTGTGTGCTG
C1	${\tt CGAGTCCCATGGCGCCGGATCGACG}$
C2	GCGTCAGGCAGC

• Building block 0 section

smiles	codon	reaction	reactant	product
OC(=O)C1=CC(=CN=C1)C#C	GCCTCG	CuAAC	$scaffold_1$	product_1
BrC1=NC=C(OCC#C)C=C1	TCCGAC	CuAAC	$scaffold_1$	$product_1$
CNC1=CC=C(OCC#C)C=C1	CAAGTG	CuAAC	$scaffold_1$	$product_1$
NC(=O)C1=CC(=CN=C1)C#C	GTCCGC	CuAAC	$scaffold_1$	$product_1$
O=C(NCC#C)NC1CC1	GACGAC	CuAAC	scaffold_1	$product_1$

\bullet Building block 1 section

smiles	codon	reaction	reactant	product
OB(O)c1cc(ccc1Cl)C#N	GTCTCAC	Suz	product_1	product_2
COc1ccc(B(O)O)c(F)c1	GTCGTAC	Suz	$product_1$	product_2

smiles	codon	reaction	reactant	product
O.Nc1cccc(c1)B(O)O	CTCATTG	Suz	$product_1$	$product_2$
Cc1cc(ccc1F)B(O)O	GTAGAGA	Suz	$product_1$	$product_2$
COc1ccc(cc1)B(O)O	GTTACCT	Suz	$product_1$	$product_2$

To facilitate the initialization of the configuration files, this information can be read in from an excel file with sheets that have the corresponding names and then converted to a config.yaml with the command

```
delt-cli init --excel_path=path/to/library.xlsx`
```

2. Enumerate library compounds:

```
Timing: 5 minutes - 2h
```

```
delt-cli library enumerate --config_path=config.yaml
```

This will create the following files in the save_dir / name folder:

- library.parquet: contains a column SMILES and the codon indices for each of the compounds.
- reaction_graph.png: a graphical representation of the reaction steps defined in the configuration file.
- 4. Calculate molecular properties: Computes comprehensive molecular descriptors including drug-likeness metrics, physical properties, and structural features.

```
delt-cli library properties --config_path=config.yaml
```

- 3. Define demultiplexing parameters: Configure the structure section of the config.yaml to enable demultiplexing with region specific error rates or allowing indels. Refer to the cutadapt documentation to learn more about calibration of the error rate.
- 4. Generate molecular representations (optional)

```
# Morgan fingerprints for similarity analysis
delt-cli library represent --method=morgan --config_path=config.yaml
# BERT embeddings for machine learning applications
delt-cli library represent --method=bert --config_path=config.yaml
```

1.9.1 Sequence demultiplexing and processing • TIMING 30 min - 4 hours

10. Execute demultiplexing workflow

```
delt-cli demultiplex run --config_path=config.yaml
```




Figure 1: reaction_graph.png




Figure 2: prop_mw.png

This performs sequential adapter trimming and barcode identification using optimized Cutadapt workflows.

11. Generate quality control reports

```
delt-cli demultiplex report --config_path=config.yaml
delt-cli demultiplex qc --config_path=config.yaml
```

Cutadapt Pipeline Report							
Region	Input	With adapters	Discarded	% with	% discarded		
0-S0	10,000	8,750	1,250	87.50%	12.50%		
1-C0	8,750	867	7,883	9.91%	90.09%		
2-B0	867	748	119	86.27%	13.73%		
3-C1	748	522	226	69.79%	30.21%		
9 -B1	522	482	40	92.34%	7.66%		
5-C2	482	448	34	92.95%	7.05%		
6-81	448	441	7	98.44%	1.56%		
Overall:							
With adapters : 441 (4.41%) Discarded : 9.559 (95.59%)							

 $hits_0-S0.pdf$

1.9.2 Statistical analysis and hit identification • TIMING 10-30 min

12. Define analysis groups

13. Perform enrichment analysis

1.9.3 Data visualization and interpretation

TIMING: 15-60 min

14. Launch interactive dashboard

Opens a web-based interface for interactive data exploration, hit visualization, and property analysis.

1.10 Troubleshooting

1.10.1 Expected performance metrics

Successful analysis should show:

- Demultiplexing efficiency >70% for high-quality data
- Even distribution of reads across expected barcode combinations
- Clear enrichment patterns in target vs control comparisons
- Consistent results across biological replicates

1.11 Timing

Protocol execution time depends on dataset size and computational resources:

- Environment setup: 15-30 minutes (one-time)
- Project initialization: 5-10 minutes
- Library enumeration: 10-60 minutes (depending on library size)
- Sequence demultiplexing: 30 minutes 4 hours (depending on read count)
- Statistical analysis: 10-30 minutes
- Visualization and interpretation: 15-45 minutes

Total workflow time: 1-6 hours for typical datasets

1.12 Anticipated results

1.12.1 Output structure

DELT-Core generates a comprehensive, standardized output structure:

```
my_decl_project/
config.yaml
```

1.12.2 Chemical library outputs

The library enumeration generates:

- Complete structure catalog: All possible compounds with canonical SMILES
- Molecular properties: Comprehensive descriptor tables for drug-likeness assessment
- Reaction network: Visual representation of synthetic pathways
- Quality validation: Structure validation and diversity analysis

For a typical 2-cycle library with 1000 building blocks per position, expect ~1 million unique structures with associated metadata.

1.12.3 Demultiplexing results

Successful processing produces:

- Barcode count tables: Quantified reads for each library member across selections
- Quality metrics: Error rates, recovery efficiency, and coverage statistics
- Statistical validation: Replicate consistency and batch effect assessment

Performance expectations:

- Initial adapter matching: 80-95% of input reads retained
- $\bullet\,$ Final barcode assignment: 60-85% of reads with valid combinations
- Library coverage: 10-90% of theoretical compounds detected
- Error rates: <5\% per position for high-quality sequencing

1.12.4 Statistical analysis outputs

The enrichment analysis provides:

- **Hit ranking**: Statistical significance and fold-change metrics
- Visualization: Volcano plots and correlation analysis
- Quality control: Replicate consistency and batch effect assessment
- Export formats: CSV tables compatible with downstream tools

1.12.5 Integration capabilities

DELT-Core outputs integrate seamlessly with:

- Machine learning workflows: Standardized feature matrices and molecular representations
- Cheminformatics pipelines: RDKit-compatible structure formats
- Statistical software: R and Python pandas-compatible data tables
- Visualization tools: Interactive dashboards and publication-ready plots

1.13 Data availability

Software, documentation, and example datasets are freely available:

- **DELT-Core repository**: https://github.com/DELTechnology/delt-core
- Documentation: Comprehensive guides and tutorials in repository wiki
- Example datasets: Test data for workflow validation
- Configuration templates: Pre-configured files for common use cases

1.14 Code availability

DELT-Core is released under the MIT License, enabling free use and modification. Complete source code with documentation is available on GitHub.

1.15 Acknowledgements

We thank the DECL research community for valuable feedback during development. We acknowledge the developers of Cutadapt, RDKit, edgeR, and other open-source tools that enable this framework.

1.16 References

[To be completed with specific citations]

1.17 Author contributions

[To be completed based on actual contributions]

1.18 Competing interests

The authors declare no competing financial interests.